

OptiGlide™ Injector



STERILE R

SINGLE USE ONLY

ISO 13485

THE USIOL ADVANTAGE

State of the art molding techniques provide precision IOL injection through 2.2 - 2.8mm incision

Fully compatible with both hydrophobic and hydrophilic acrylic lenses

Designed to accommodate thin lens technology for micro incision technique

Side windows allow for surgeon's inspection before implantation

Provides smooth, controlled passage through cartridge, eliminating IOL scratches, haptic stretching and incision widening

| | |
|----------------|---------------|
| Tip shape | Bevel down |
| Outer diameter | 1.6mm x 2.0mm |
| Angulation | 25° |

S H A R I N G Y O U R V I S I O N

